

Standard Analysis Certificate

PRODUCT :n-HEPTANE 95% RS HPLC

CODE :446830

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20° C	-	0.681 - 0.687
Refractive index at 20°C	-	1.386 - 1.390
Water (K.F.)	ppm	<= 100
Residue on evaporation	ppm	<= 5
Acidity or alkalinity	meq/g	<= 0.00015
Assay (GLC)	%	>= 95
U.V. Transmittance	-	-
at 200 nm	%	>= 20
at 210 nm	%	>= 45
at 220 nm	%	>= 80
At 250 nm	%	>= 98
Aromatic compounds	ppm	<= 10

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)