

# Standard Analysis Certificate

**PRODUCT** :ISOHEXANE RS PESTIPUR - For pesticide analysis  
**CODE** :447130  
**METHOD** :LAB118

| TEST                                     | U.M.  | SPECIFICATION |
|--|-------|---------------|
| Identification (IR)                      | -     | Conform       |
| Water content (K.F.)                     | mg/Kg | <= 100        |
| Non volatile residue                     | mg/Kg | <= 5          |
| Colour                                   | Hazen | <= 10         |
| Total C6H14 isomers (GC)                 | %     | >= 80         |
| n-hexane                                 | %     | <= 5          |
| GC-ECD.Individual peak (Lindane)         | ng/l  | <= 3          |
| Retention time trichlorobenzene to mirex | -     | -             |

For chlorinated compounds analysis; extraction of pesticides.

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

**QUALITY CONTROL RESPONSIBLE**

**B. COULANGE (VDR)**