

# Standard Analysis Certificate

**PRODUCT** :DIETHYLETHER RS SPECTROSOL Not stabilized  
**CODE** :447590  
**METHOD** :9903

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Colour (APHA)	-	<=10
Identification	-	Positive
Density at 20° C	-	0.712 ÷ 0.714
Boiling point	° C	34.4 ÷ 34.9
Water (K.F.)	ppm	<=100
Residue on evaporation	ppm	<=5
Acidity	meq/g	<=0.0005
Alcalinity	meq/g	<=0.0002
Peroxides (H2O2)	ppm	<=60
Assay (GLC)	%	>=99.7
U.V. Transmittance	-	-
at 220 nm	%	>=10
at 230 nm	%	>=45
at 250 nm	%	>=75
at 280 nm	%	>=92
at 300 nm	%	>=98

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)