

TECHNICAL DATA SHEET

| | |
|--------------------|--------|
| Part Number | 451000 |
| Description | Stool |



Benefits :

- ✓ TABAT sit-stand stool
- ✓ Base and seat support in 304 stainless steel, seat in molded polyurethane
- ✓ Created especially for clean rooms
- ✓ Branches in flat stainless steel and not in tubes, do not have any retention zone
- ✓ Perfect stability thanks to the 5 feet mounted on stainless steel jacks
- ✓ Longevity and robustness

Technical specifications :

- ✓ Height adjustment by stainless steel screws from 490 to 790 mm
- ✓ Base \varnothing 530 mm
- ✓ Weight 11 kg

Sales unit :

1 Stool
