

EC Declaration of Conformity

Conformity to DIRECTIVE 98/79/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 27 October 1998 on in vitro diagnostic medical devices

Manufacturer: Greiner Bio-One GmbH
Bad Haller Straße 32
4550 Kremsmünster
Austria

Production Greiner Bio-One GmbH
Location: Bad Haller Straße 32
4550 Kremsmünster
Austria

Greiner Bio-One North America Inc.
4238 Capital Drive, Monroe
NC 28110
United States of America

Greiner Bio-One Brasil
Produtos Medicos
Hospitalares Ltda.
Av. Affonso Pansan no. 1.967
13473-620, Villa Bertini
Americana, Sao Paulo
Brasil

Greiner Bio-One (Thailand) Ltd.
700/172 Moo.1
Amata Nakorn Industrial Estate
Tambon Bankao
Amphur Phanthong
Chonburi 20160
Thailand

Product / VACUETTE® TUBES
Product Group: (for details please refer to page 2-33)

Classification: Other device (all devices except Annex II and except self-testing devices)

GMDN Code(s): (for details please refer to page 2-33)

We herewith declare under our sole responsibility that the above mentioned products meet the provisions of the above EC Council Directive and the applicable standards. All supporting documentations are retained under the premises of the manufacturer.

Conformity Assessment procedure acc. to Annex III of the Directive 98/79/EC of the European Parliament and of the Council on in vitro diagnostic medical devices.

Standards:

Refer to the List of applicable (harmonized) standards in the Technical Documentation.

Kremsmünster, 26.07.2022



Signature:

Georg Sambs
Quality Manager GBO AT



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2 ml 4NC ESR sodium citrate 3.2% 13x75 black cap-black ring, PET, PREMIUM	454073	47592
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2 ml 4NC ESR sodium citrate 3.2% 13x75 black cap-black ring, G-barcode label, PREMIUM	484516	47592
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2.75 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, PP	729070	47592
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.5 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, PP	729073	47592
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.5 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, PP, High Altitude	729073H	47592
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2.9 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, glass	729090	47592
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.6 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, glass	729093	47592
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.6 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, G-barcode label, glass	769093	47592
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.5 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, PP, G-barcode label	769094	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, PREMIUM	454321	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged	454322	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged, High Altitude	454322H	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged, China	454322CN	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged, Korea	454322KR	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged, High Altitude	454323	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, PREMIUM	454325	47592



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, PREMIUM	454327	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 4.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, PREMIUM	454329	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 light blue cap-black ring, sandwich tube, non-ridged, High Altitude	454331	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, non-ridged	454332	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, special label, non-ridged, CSS	454332CSS	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, transparent label, non-ridged	454333	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, transparent special label, non-ridged	454333CN	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, non-ridged	454334	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3ml 9NC Coagulation sodium citrate 3.2% 13x75 light blue cap-black Ring, sandwich tube, non-ridged, Canada	454334CA	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, non-ridged, China	454334CN	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, non-ridged, High-Altitude	454337H	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3.2% 13x75 light blue cap-black ring, sandwich tube, non-ridged, High Altitude	454351	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged	454352	47592

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, transparent label, non-ridged	454360	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 1.8 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, transparent special label, non-ridged, Japan	454360JP	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, transparent label, non-ridged	454471	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 9 ml 9NC Coagulation sodium citrate 3,2% 16x100 blue cap-black ring, non-ridged	455322	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, transparent label, PREMIUM	474321	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 1.8 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, transparent special label, PREMIUM, Japan	474321JP	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 1.8 ml 9NC Coagulation sodium citrate 3.2% 13x75 black cap-white ring, sandwich tube, transparent special label, PREMIUM, Japan	474322JP	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, transparent label, PREMIUM	474327	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, ISLAB barcode label, sandwich tube, non-ridged	484322	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, ISLAB barcode label M, non-ridged	484322M	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, KOH barcode label, non-ridged	484334	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, G-barcode label, non-ridged	484508	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, G-barcode label, PREMIUM	484517	47592

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, G-barcode label, PREMIUM	484518	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 5 ml 9NC Coagulation sodium citrate 3.2% 13x100 blue cap-black ring, G-barcode label, PREMIUM	486505	47592
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, G-bi-barcode label, sandwich tube, PREMIUM	494512	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-white ring, sandwich tube, PREMIUM	454381	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-white ring, sandwich tube, PREMIUM, China	454381CN	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-white ring, sandwich tube, special label, PREMIUM, High Altitude	454383CN	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-white ring, sandwich tube, special label, PREMIUM, High Altitude	454383CNH	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, sandwich tube, PREMIUM	454385	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, sandwich tube, PREMIUM, China	454385CN	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, sandwich tube, PREMIUM	454387	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 4.5 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, PREMIUM	454389	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, sandwich tube, non-ridged	454392	47592
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 9 ml 9NC Coagulation sodium citrate 3,8% 16x100 blue cap-black ring, non-ridged	455382	47592
CAT	Crossmatch Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Crossmatch Serum Clot Activator 13x100 pink cap-black ring, special label non-ridged	456097	42386



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Clot Activator	VACUETTE® TUBE 4.5 ml CAT Serum Clot Activator 13x75 red cap-black ring, PREMIUM	454027	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, PREMIUM	454092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 3 ml CAT Serum Clot Activator 13x75 red cap-black ring, PREMIUM	454095	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 2 ml CAT Serum Clot Activator 13x75 red cap-white ring, PREMIUM	454096	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 1 ml CAT Serum Clot Activator 13x75 red cap-white ring, PREMIUM	454098	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 white cap-black ring, PREMIUM	454200	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 pink cap-black ring, PREMIUM	454202	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, non-ridged	454204	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, non-ridged, High Altitude	454204H	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, non-ridged, China	454204CN	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 2 ml CAT Serum Clot Activator 13x75 red cap-white ring, non-ridged	454236	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 pink cap-black ring, non-ridged	454262	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, transparent label, non-ridged	454433	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 2 ml CAT Serum Clot Activator 13x75 red cap-white ring, transparent label, non-ridged	454461	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, transparent label, non-ridged	455009	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, non-ridged	455092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, non-ridged, High Altitude	455092H	42386

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, special label, non-ridged, Japan	455092JP	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, special label, non-ridged	455092CSS	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 gold cap-gold ring, non-ridged	455204	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 pink cap-black ring, pink label, PREMIUM	456014	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 gold cap-black ring, black label, PREMIUM	456015	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, transparent label, non-ridged	456048	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, non-ridged	456089	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, non-ridged, High Altitude	456089H	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, non-ridged, Korea	456089KR	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, special label, non-ridged, High Altitude	456089CNH	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 orange cap-black ring, non-ridged	456091	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, PREMIUM	456092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 pink cap-black ring, non-ridged	456248	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-red ring, non-ridged	456249	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4.5 ml CAT Serum Clot Activator 13x75 red cap-black ring, transparent label, PREMIUM	474027	42386



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, transparent label, PREMIUM	474092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, transparent label, PREMIUM	476092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, transparent special label, PREMIUM, Japan	476092JP	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 2 ml CAT Serum Clot Activator 13x75 red cap-white ring, G-barcode label, PREMIUM	484515	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, G-barcode label, PREMIUM	484525	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, G-barcode label, non-ridged	485504	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, ISLAB barcode label, non-ridged	486089	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, ISLAB barcode label M, non-ridged	486089M	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, G-barcode label, PREMIUM	486504	42386
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM	454028	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM, China	454028CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged	454067	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged, USA	454067P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged	454067R	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged, China	454067CN	41128



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged, High Altitude, USA	454067PH	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM	454071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM, China	454071CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM, Hong Kong	454071HK	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, special label, PREMIUM, OMB	454071OMB	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, transparent label, non-ridged	454078	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, transparent label, non-ridged, P	454078P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, special label, non-ridged	454078CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 yellow cap-white ring, black label, non-ridged	454205	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 pink cap-yellow ring, pink label, PREMIUM	454212	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, PREMIUM	454214	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, transparent special label, non-ridged	454219CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged	454225	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, black label, non-ridged	454228	41128



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, black label, non-ridged, USA	454228P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, black label, non-ridged, China	454228CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, non-ridged	454243	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, non-ridged, USA	454243P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, non-ridged, China	454243CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-white ring, black label, non-ridged	454253	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 white cap-yellow ring, PREMIUM	454260	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 royal blue cap-yellow ring, non-ridged	454263	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 yellow cap-yellow ring, non-ridged	454264	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-black ring, non-ridged	454319	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, transparent label, non-ridged	454420	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, transparent label, non-ridged, China	454420CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, transparent label, non-ridged	454473	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, photoresistant, non-ridged	454482	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455010	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, USA	455010P	41128



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 gold cap-gold ring, non-ridged	455034	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455057	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, Gamma	455071G	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455071P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455071R	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, special label, non-ridged, Korea	455071KR	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 gold cap-gold ring, non-ridged	455072	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 white cap-black ring, black label, non-ridged	455078	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 7 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, double gel	455099	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 7 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, double gel, USA	455099P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, transparent label, non-ridged	455211	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, transparent special label, non-ridged, Korea	455211KR	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 white cap-yellow ring, black label, non-ridged	456004	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-yellow ring, black label, non-ridged	456005	41128



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-yellow ring, black special label, non-ridged	456005CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 royal blue cap-yellow ring, blue label, non-ridged	456006	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, PREMIUM	456010	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black label, non-ridged	456018	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black label, non-ridged	456018P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black special label, non-ridged	456018CA	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black label, non-ridged, China	456018CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black label, non-ridged, High Altitude, USA	456018PH	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black special label, non-ridged	456018PCA	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 pink cap-yellow ring, non-ridged	456019	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent label, non-ridged	456069	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent special label, non-ridged	456069CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent special label, non-ridged	456069JP	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent label, non-ridged, Korea	456069KR	41128



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, PREMIUM	456071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, PREMIUM	456071P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, PREMIUM, China	456071CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, special label, PREMIUM, High Altitude	456071CNH	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged	456073	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged	456073P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged	456073R	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged, China	456073CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged, Korea	456073KR	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged, High Altitude, USA	456073PH	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, special label, non-ridged	456073PCSS	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-black ring, black label, PREMIUM	456077	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 white cap-black ring, black label, PREMIUM	456078	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 orange cap-black ring, black label, PREMIUM	456079	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 royal blue cap-yellow ring, transparent special label, PREMIUM, Japan	456212JP	41128



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, transparent label, non-ridged	456234	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, transparent label, non-ridged, China	456234CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-red ring, non-ridged	456239	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 white cap-white ring, black label, non-ridged	456278	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 pink cap-black ring, special label, non-ridged	456288	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-black ring, special label, non-ridged	456289	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-black ring, special label, non-ridged	456290	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged, USA	456292P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator Lipo Tube 13x100 black cap-yellow ring, non-ridged, USA	456293P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-white ring, non-ridged	456295	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, transparent label, PREMIUM	474071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, transparent label, PREMIUM	474214	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, transparent label, PREMIUM	476010	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, transparent special label, with circle fill line, PREMIUM, Japan	476010JP	41128



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent label, PREMIUM	476071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent special label, PREMIUM, Japan	476071JP	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 yellow cap-white ring, black KOH barcode label, non-ridged	484485	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, G-barcode label, non-ridged	484502	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, G-barcode label, PREMIUM	484513	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, G-barcode label, PREMIUM	484514	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, barcode label, non-ridged	485071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator 16x100 red cap-yellow ring, G-barcode label, non-ridged	485502	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, G-barcode label, non-ridged	485509	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, G-barcode label, non-ridged	485514	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 white cap-yellow ring, KOH barcode label, non-ridged	486004	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, KOH barcode label, non-ridged	486005	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 royal blue cap-yellow ring, KOH barcode label, non-ridged	486006	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-yellow ring, KOH barcode label, non-ridged	486007	41128



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black KOH barcode label, non-ridged	486018	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 pink cap-yellow ring, KOH barcode label, non-ridged	486019	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, ISLAB barcode label, non-ridged	486073	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, ISLAB barcode label M, non-ridged	486073M	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black G-barcode label, non-ridged	486500	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, G-barcode label, PREMIUM	486503	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, G-bi-barcode label, non-ridged	495071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, G-bi-barcode label, PREMIUM	496503	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 3.5 ml CAT Serum Fast Separator 13x75 orange cap-yellow ring, non-ridged	454591	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 3.5 ml CAT Serum Fast Separator 13x75 orange cap-yellow ring, PREMIUM	454592	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 3.5 ml CAT Serum Fast Separator 13x75 orange cap-yellow ring, non-ridged	454593	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 5 ml CAT Serum Fast Separator 13x100 orange cap-yellow ring, PREMIUM	456309	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 5 ml CAT Serum Fast Separator 13x100 orange cap-yellow ring, special label, PREMIUM, Japan	456309JP	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 5 ml CAT Serum Fast Separator 13x100 orange cap-yellow ring, non-ridged	456310	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 5 ml CAT Serum Fast Separator 13x100 orange cap-yellow ring, non-ridged	456313	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 5 ml CAT Serum Fast Separator 13x100 orange cap-yellow ring, special label, PREMIUM, Japan	476309JP	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 5 ml CAT Serum Fast Separator 13x100 orange cap-yellow ring, non-ridged	486509	41128
FE	Glucose Test only	VACUETTE® TUBE 2 ml FE Glucose Test only 13x75 grey cap-white ring, non-ridged	454252	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 4 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-black ring, non-ridged	454033	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, PREMIUM	454085	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 4 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-black ring, PREMIUM	454091	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, non-ridged	454221	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, transparent label, PREMIUM	454344	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 4 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-black ring, KOH barcode label, non-ridged	484033	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, ISLAB barcode label, non-ridged	484221	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, ISLAB barcode label M, non-ridged	484221M	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, G-barcode label, PREMIUM	484520	42329
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 4 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-black ring, G-barcode label, PREMIUM	484527	42329
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, PREMIUM	454061	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, PREMIUM, Japan	454061JP	47591



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, PREMIUM	454062	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, transparent label, non-ridged	454231	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, non-ridged	454238	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, non-ridged	454297	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-yellow ring, non-ridged	454427	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, transparent label, non-ridged	454438	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 5,5 ml FX Sodium Fluoride / Potassium Oxalate 13x100 grey cap-black ring, non-ridged	456020	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 6 ml FX Sodium Fluoride / Potassium Oxalate 13x100 grey cap-black ring, PREMIUM	456062	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, G-barcode label, PREMIUM	484519	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, G-barcode label, non-ridged	484528	47591
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged	454009	47592
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, PREMIUM	454020	47592
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, special label, PREMIUM, High Altitude	454020CNH	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, PREMIUM	454023	47592
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, PREMIUM	454024	47592
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, PREMIUM, China	454024CN	47592



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, special label, PREMIUM, Japan	454024JP	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-yellow ring, PREMIUM	454032	47592
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, non-ridged	454047	47592
K2E	K2EDTA	VACUETTE® TUBE 1 ml K2E K2EDTA 13x75 lavender cap-white ring, non-ridged	454052	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 pink cap-black ring, non-ridged	454208	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged	454209	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged, High Altitude	454209H	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged, Canada	454209CA	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged, China	454209CN	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-white ring, PREMIUM	454210	47592
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-red ring, non-ridged	454234	47592
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged	454246	47592
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, special label, non-ridged	454246CA	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-white ring, non-ridged	454286	47592
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 purple cap-black ring, non-ridged	454346	47592
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, transparent label, non-ridged	454409	47592
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, transparent special label, non-ridged, Japan	454409JP	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, transparent label, non-ridged	454410	47592
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, transparent label, non-ridged	454411	47592
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, non-ridged, USA	454428	47592
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 yellow cap-black ring, non-ridged	454504	47592
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 yellow cap-white ring, non-ridged	454506	47592



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, barcode label, non-ridged	484246	47592
K2E	K2EDTA	VACUETTE® TUBE 9 ml K2E K2EDTA 16x100 lavender cap-black ring, non-ridged	455045	47592
K2E	K2EDTA	VACUETTE® TUBE 9 ml K2E K2EDTA 16x100 black cap-black ring, non-ridged	455075	47592
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 lavender cap-black ring, non-ridged, USA	456002	47592
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 lavender cap-black ring, PREMIUM	456023	47592
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 lavender cap-black ring, non-ridged	456043	47592
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 pink cap-black ring, transparent special label, non-ridged	456068	47592
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 lavender cap-black ring, transparent label, non-ridged	456243	47592
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 red cap-black ring, PREMIUM	456244	47592
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 pink cap-pink ring, non-ridged	456279	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 pink cap-black ring, G-barcode label, non-ridged	484506	47592
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, G-barcode label, non-ridged	484507	47592
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, G-barcode label, PREMIUM	484526	47592
K2E	K2EDTA	VACUETTE® TUBE 9 ml K2E K2EDTA 16x100 lavender cap-black ring, G-Bi-Barcode label, non-ridged	495045	47592
K2E	K2EDTA Sep	VACUETTE® TUBE 3.5 ml K2E K2EDTA Separator 13x75 lavender cap-yellow ring, PREMIUM	454235	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 3.5 ml K2E K2EDTA Separator 13x75 lavender cap-yellow ring, PREMIUM	454235P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 8 ml K2E K2EDTA Separator 16x100 white cap-yellow ring, non-ridged	455021	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 8 ml K2E K2EDTA Separator 16x100 lavender cap-yellow ring, non-ridged	455040	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 8 ml K2E K2EDTA Separator 16x100 lavender cap-yellow ring, non-ridged, USA	455040P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 8 ml K2E K2EDTA Separator 16x100 lavender cap-yellow ring, non-ridged, China	455040CN	47587



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, PREMIUM	456011	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, PREMIUM, China	456011CN	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 4.5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, non-ridged, USA	456024P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, non-ridged	456058	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, non-ridged	456058P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, special label, non-ridged	456058CN	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, transparent label, non-ridged	456205P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, PREMIUM	456211	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, special label, non-ridged	456225	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, non-ridged	456282	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, non-ridged	456282P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, special label, non-ridged, Japan	456282JP	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, G-barcode label, Premium	486511P	47587
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, non-ridged	456052	47592
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, PREMIUM	456061	47592



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, PREMIUM	456093	47592
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, non- ridged	456152	47592
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, non- ridged	456242	47592
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, non- ridged	456252	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-red ring, PREMIUM	454003	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 pink cap-black ring, non-ridged	454012	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-red ring, non-ridged	454015	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged	454021	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged, High Altitude	454021H	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 red cap-black ring, non-ridged	454022	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 orange cap-black ring, non-ridged	454025	47592
K3E	K3EDTA	VACUETTE® TUBE 1 ml K3E K3EDTA 13x75 lavender cap-white ring, PREMIUM	454034	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, PREMIUM	454036	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, PREMIUM, China	454036CN	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA Liquid 13x75 lavender cap-black ring, PREMIUM	454041	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, PREMIUM	454086	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, PREMIUM	454087	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, PREMIUM, China	454087CN	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 ping cap-black ring, pink label, PREMIUM	454099	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged	454217	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged, High Altitude	454217H	47592



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged, Korea	454217KR	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, non-ridged	454222	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, non-ridged, High Altitude	454222H	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, non-ridged, China	454222CN	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-white ring, PREMIUM	454232	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 black cap-black ring, PREMIUM	454245	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, transparent label, non-ridged	454287	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, transparent label, non-ridged, Korea	454287KR	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent label, non-ridged	454317	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent special label, non-ridged	454317CN	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged 100/Box	454408	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-red ring, transparent label, non-ridged	454515	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, amber tube, non-ridged	454921	47592
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 black cap-black ring, black label, non- ridged	455030	47592
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 lavender cap-black ring, non-ridged	455036	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-black ring, non-ridged	456003	47592
K3E	K3EDTA	VACUETTE® TUBE 5 ml K3E K3EDTA 13x100 lavender cap-black ring, PREMIUM	456034	47592
K3E	K3EDTA	VACUETTE® TUBE 5 ml K3E K3EDTA 13x100 lavender cap-black ring, PREMIUM, China	456034CN	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, PREMIUM	456036	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA Liquid 13x100 lavender cap-black ring, PREMIUM	456037	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, non-ridged	456038	47592



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, non-ridged, Taiwan	456038TW	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 yellow cap-black ring, PREMIUM	456039	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 royal blue cap-black ring, non-ridged	456042	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, transparent label, non-ridged	456053	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-black ring, transparent label, non-ridged	456067	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 white cap-black ring, non-ridged	456074	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 white cap-black ring, non-ridged, Taiwan	456074TW	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-black ring, pink label, PREMIUM	456099	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 green cap-black ring, non-ridged	456236	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-red ring, special label, non-ridged	456251	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-red ring, special label, non- ridged	456257	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-yellow ring, special label, non- ridged	456259	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 yellow cap-red ring, special label, non- ridged	456260	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 yellow cap-yellow ring, special label, non-ridged	456262	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 red cap-red ring, special label, non- ridged	456263	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 red cap-yellow ring, special label, non- ridged	456266	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, transparent label, PREMIUM	474036	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, transparent label, PREMIUM	474086	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent label, PREMIUM	474087	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent label, PREMIUM	474232	47592



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 orange cap-black ring, KOH barcode label, non-ridged	484025	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, KOH barcode label, non-ridged	484217	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, ISLAB barcode label, non-ridged	484317	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, ISLAB barcode label M, non-ridged	484317M	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, KOH barcode label, non-ridged	484481	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, G-barcode label, non-ridged	484500	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, G-barcode label, non-ridged	484501	47592
K3E	K3EDTA	VACUETTE® TUBE 1 ml K3E K3EDTA 13x75 lavender cap-white ring, G-barcode label, PREMIUM	484509	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, G-barcode label, PREMIUM	484510	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, G-barcode label, PREMIUM	484511	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, G-barcode label, PREMIUM	484512	47592
K3E	K3EDTA	VACUETTE® TUBE 4.5 ml K3E K3EDTA 13x75 lavender cap-black ring, G-barcode label, non-ridged	484529	47592
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 lavender cap-black ring, G-barcode label, non-ridged	485501	47592
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 orange cap-black ring, G-barcode label, non-ridged	485512	47592
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, G-barcode label, PREMIUM	486502	47592
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, G-bi-barcode label, non-ridged	494217	47592
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, G-bi-barcode label, non-ridged	494508	47592
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 ping cap-black ring, pink G-bi-barcode label, PREMIUM	494509	47592



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, G-bi-barcode label, non-ridged	496511	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged	454029	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged, High Altitude	454029H	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged, China	454029CN	47592
LH	Lithium Heparin	VACUETTE® TUBE 4.5 ml LH Lithium Heparin 13x75 green cap-black ring, PREMIUM	454049	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, transparent label, non-ridged	454056	47592
LH	Lithium Heparin	VACUETTE® TUBE 1 ml LH Lithium Heparin 13x75 green cap-white ring, PREMIUM	454081	47592
LH	Lithium Heparin	VACUETTE® TUBE 3 ml LH Lithium Heparin 13x75 green cap-black ring, PREMIUM	454082	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, PREMIUM	454084	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, special label, PREMIUM	454084CN	47592
LH	Lithium Heparin	VACUETTE® TUBE 2 ml LH Lithium Heparin 13x75 green cap-white ring, PREMIUM	454089	47592
LH	Lithium Heparin	VACUETTE® TUBE 2 ml LH Lithium Heparin 13x75 green cap-white ring, PREMIUM, China	454089CN	47592
LH	Lithium Heparin	VACUETTE® TUBE 2 ml LH Lithium Heparin 13x75 green cap-white ring, non-ridged	454237	47592
LH	Lithium Heparin	VACUETTE® TUBE 3 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged	454244	47592
LH	Lithium Heparin	VACUETTE® TUBE 3 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged, China	454244CN	47592
LH	Lithium Heparin	VACUETTE® TUBE 9 ml LH Lithium Heparin 16x100 green cap-black ring, non-ridged	455084	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 16x100 green cap-black ring, non-ridged	455222	47592
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, PREMIUM	456084	47592
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, non-ridged	456088	47592
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, special label, non-ridged, China	456088CN	47592
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, transparent label, non-ridged	456271	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, transparent label, PREMIUM	474084	47592

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, transparent label, PREMIUM	476084	47592
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, transparent special label, PREMIUM, Japan	476084JP	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, KOH barcode label, non-ridged	484029	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, G-barcode label, PREMIUM	484522	47592
LH	Lithium Heparin	VACUETTE® TUBE 2 ml LH Lithium Heparin 13x75 green cap-white ring, G-barcode label, PREMIUM	484589	47592
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 16x100 green cap-black ring, G-barcode label, non-ridged	485222	47592
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, G-barcode label, PREMIUM	486508	47592
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged	454008	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged	454008P	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged, China	454008CN	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 2.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, PREMIUM	454046	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, PREMIUM	454083	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged	454247	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged	454247P	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 gold cap-white ring, non-ridged	454446	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 gold cap-white ring, non-ridged, USA	454446P	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, special label, non-ridged	454483	47587



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 mint green cap-yellow ring, transparent label, non-ridged	454489	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 mint green cap-yellow ring, non-ridged, special label	454490PCA	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, transparent label, non-ridged	454491P	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 mint green cap-yellow ring, transparent label, non-ridged	454493	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 mint green cap-yellow ring, non-ridged	454493P	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 8 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, non-ridged	455083	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 8 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, non-ridged	455083P	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 4 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, non-ridged	455229	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, PREMIUM	456083	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, PREMIUM	456083P	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, PREMIUM, China	456083CN	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged	456087	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged	456087P	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged, Korea	456087KR	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, transparent label, non-ridged	456216	47587



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, transparent special label, non-ridged, Japan	456216JP	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, transparent label, non-ridged, Korea	456216KR	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged	456287	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged, USA	456287P	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 mint green cap-yellow ring, special label, non-ridged	456289CA	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 mint green cap-yellow ring, special label, non-ridged	456289PCA	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-white ring, non-ridged	456303	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 mint green cap-yellow ring, non-ridged	456305	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, transparent label, PREMIUM	474080	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, transparent label, PREMIUM	474083	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, transparent label, PREMIUM	476083	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-red ring, barcode label, non-ridged	484030	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, barcode label, non-ridged	484031	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, G-barcode label, non-ridged	484523	47587



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
LH	Lithium Heparin Sep	VACUETTE® TUBE 4 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, G-barcode label, non-ridged	485513	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, ISLAB barcode label, non-ridged	486087	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, ISLAB barcode label M, non-ridged	486087M	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, G-barcode label, non-ridged	486507	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, G-bi-barcode label, non-ridged	494504	47587
LH	Lithium Heparin Sep	VACUETTE® TUBE 8 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, G-Bi-Barcode label, non-ridged	495083	47587
LH	Lithium Heparin/ Iodoacetate	VACUETTE® TUBE 4.5 ml Lithium Heparin / Iodoacetate 13x75 grey cap-black ring, PREMIUM	454080	47592
LH	Lithium Heparin/ Iodoacetate	VACUETTE® TUBE 4.5 ml Lithium Heparin / Iodoacetate 13x75 grey cap-black ring, non-ridged	454248	47592
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 green cap-green ring, non-ridged	454030	47592
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 green cap-green ring, non-ridged, High Altitude	454030H	47592
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 green cap-green ring, PREMIUM	454051	47592
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 green cap-green ring, special label, PREMIUM	454051CN	47592
NH	Sodium Heparin	VACUETTE® TUBE 2 ml NH Sodium Heparin 13x75 green cap-green ring, non-ridged	454302	47592
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 royal blue cap-black ring, transparent label, non-ridged	454502	47592
NH	Sodium Heparin	VACUETTE® TUBE 9 ml NH Sodium Heparin 16x100 green cap-green ring, non-ridged	455051	47592
NH	Sodium Heparin	VACUETTE® TUBE 6 ml NH Sodium Heparin 13x100 green cap-green ring, non-ridged	456028	47592
NH	Sodium Heparin	VACUETTE® TUBE 6 ml NH Sodium Heparin 13x100 green cap-green ring, special label, non-ridged, Japan	456028JP	47592
NH	Sodium Heparin	VACUETTE® TUBE 6 ml NH Sodium Heparin 13x100 green cap-green ring, PREMIUM	456051	47592



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 green cap-green ring, G-barcode label,	484531	47592
NH	Sodium Heparin	VACUETTE® TUBE 9 ml NH Sodium Heparin 16x100 green cap-black ring, G-barcode label, non-ridged	485505	47592
NH	Trace Elements Sodium Heparin	VACUETTE® TUBE 6ml NH Trace Elements Sodium Heparin 13x100 royal blue cap-black ring, PREMIUM	456080	47592
NH	Trace Elements Sodium Heparin	VACUETTE® TUBE 6 ml NH Trace Elements Sodium Heparin 13x100 royalblue cap-black ring, non-ridged	456275	47592
NH	Trace Elements Sodium Heparin	VACUETTE® TUBE 6ml NH Trace Elements Sodium Heparin 13x100 royal blue cap-black ring, ISLAB barcode label, PREMIUM	486080	47592
NH	Trace Elements Sodium Heparin	VACUETTE® TUBE 6ml NH Trace Elements Sodium Heparin 13x100 royal blue cap-black ring, ISLAB barcode label M, PREMIUM	486080M	47592
Others	ACD-A	VACUETTE® TUBE 9 ml ACD-A 16x100 yellow cap-black ring, non-ridged	455055	47592
Others	ACD-A	VACUETTE® TUBE 6 ml ACD-A 13x100 yellow cap-black ring, PREMIUM	456055	47592
Others	ACD-B	VACUETTE® TUBE 9 ml ACD-B 16x100 yellow cap-black ring, non-ridged	455094	47592
Others	ACD-B	VACUETTE® TUBE 6 ml ACD-B 13x100 yellow cap-black ring, PREMIUM	456094	47592
Others	CPDA	VACUETTE® TUBE 9 ml CPDA 16x100 yellow cap-black ring, non-ridged	455056	47592
Others	CPDA	VACUETTE® TUBE 6 ml CPDA 13x100 yellow cap-black ring, PREMIUM	456056	47592
Others	CPDA	VACUETTE® TUBE 6 ml CPDA 13x100 yellow cap-black ring, non-ridged	456057	47592
Others	CTAD	VACUETTE® TUBE 3.5 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, PREMIUM	454064	47592
Others	CTAD	VACUETTE® TUBE 3.5 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, transparent label, non-ridged	454462	47592
Others	CTAD	VACUETTE® TUBE 3.5 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, transparent label, PREMIUM	474064	47592
Others	CTAD	VACUETTE® TUBE 2 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, transparent label, PREMIUM	474065	47592
Others	CTAD	VACUETTE® TUBE 3 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, transparent label, PREMIUM	474304	47592
Others	EDTA Aprotinin	VACUETTE® TUBE 4 ml EDTA Aprotinin 13x75 pink cap-black ring, PREMIUM	454261	58135
Others	FC-Mix	VACUETTE® FC Mix TUBE 2 ml 13x75 pink cap-white ring, transparent label, non-ridged	454510	47591



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
Others	FC-Mix	VACUETTE® FC Mix TUBE 2 ml 13x75 grey cap-white ring, transparent label, non-ridged	454511	47591
Others	FC-Mix	VACUETTE® FC Mix TUBE 3 ml 13x75 grey cap-black ring, non-ridged	454512	47591
Others	FC-Mix	VACUETTE® FC Mix TUBE 3 ml 13x75 pink cap-black ring, transparent label, non-ridged	454513	47591
Others	FC-Mix	VACUETTE® FC Mix TUBE 3 ml 13x75 grey cap-black ring, transparent label, non-ridged	454514	47591
Others	Homocysteine Detection Tube	VACUETTE® TUBE 2 ml Homocysteine Detection 13x75 white cap-red ring, sandwich tube, PREMIUM	454421	47592
Z	Discard Tube	VACUETTE® TUBE 5 ml Z Discard Tube 13x100 black cap-red ring, non-ridged	456202	47590
Z	Discard Tube	VACUETTE® TUBE 5 ml Z Discard Tube 13x100 white cap-black ring, non-ridged	456217	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z No Additive 13x75 white cap-black ring, PREMIUM	454001	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z No Additive 13x75 white cap-black ring, PREMIUM China	454001CN	47590
Z	No Additive	VACUETTE® TUBE 3.5 ml Z No Additive 13x75 white cap-black ring, non-ridged	454045	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 white cap-white ring, PREMIUM	454088	47590
Z	No Additive	VACUETTE® TUBE 3 ml Z No Additive 13x75 red cap-white ring, transparent label, non-ridged	454224	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 blue cap-black ring, non-ridged	454240	47590
Z	No Additive	VACUETTE® TUBE 3 ml Z No Additive 13x75 white cap-black ring, non-ridged	454241	47590
Z	No Additive	VACUETTE® TUBE 3 ml Z No Additive 13x75 gold cap-gold ring, red label, non-ridged	454256	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 black cap-black ring, non-ridged	454318	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 black cap-yellow ring, non-ridged, High Altitude	454318H	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 lavender cap-black ring, non-ridged	454370	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 lavender cap-yellow ring, non-ridged, High Altitude	454370H	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 pink cap-black ring, non-ridged	454371	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 pink cap-yellow ring, non-ridged, High Altitude	454371H	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 light blue cap-black ring, non-ridged	454372	47590



VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 light blue cap-yellow ring, non-ridged, High Altitude	454372H	47590
Z	No Additive	VACUETTE® TUBE 9 ml Z No Additive 16x100 white cap-black ring, non-ridged	455001	47590
Z	No Additive	VACUETTE® TUBE 6 ml Z No Additive 13x100 white cap-black ring, PREMIUM	456001	47590
Z	No Additive	VACUETTE® TUBE 6 ml Z No Additive 13x100 white cap-black ring, non-ridged	456085	47590
Z	No Additive	VACUETTE® TUBE 5 ml Z No Additive 13x100 white cap-black ring	456285	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z No Additive 13x100 lavender cap-black ring	456286	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 white cap-white ring, G-barcode label, PREMIUM	484088	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z No Additive 13x75 white cap-black ring, G-barcode label, PREMIUM	484530	47590
Z	No Additive	VACUETTE® TUBE 9 ml Z No Additive 16x100 white cap-black ring, G-barcode label, non-ridged	485506	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z No Additive 13x75 white cap-black ring, G-barcode label, PREMIUM	494510	47590
Z	No Additive	VACUETTE® TUBE 3 ml Z No Additive 13x75 gold cap-gold ring, red G-barcode label, non-ridged	494511	47590
Z	Trace Elements No Additive	VACUETTE® TUBE 6 ml Z Trace Elements No Additive 13x100 royal blue cap-black ring, non-ridged	456016	47590