

# Standard Analysis Certificate

**PRODUCT** :IODINE 0.1 N (0.05 mol/L) RPE**CODE** :456030000

TEST	U.M.	SPECIFICATION
Description Assay (potentiometry)	- N	Brown red liquid 0.0998 - 0.1002

This volumetric solution is checked by sodium thiosulfate solution, tested with potassium dichromate volumetric standard traceable to the SRM 136e of NIST (National Institute of Standards and Technology)

DOMINIQUE DUTSCHER SAS

**Date**

:11/01/2020

**QUALITY CONTROL RESPONSIBLE****B. COULANGE (VDR)**