

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	46270-50ML-HON
<b>Description</b>	Ferrouin indicator solution

---

## Technical specifications :

- ✓ ferrouin indicator solution, redox indicator, 0,025 m 50ml
- ✓ CAS N°14634-91-4
- ✓ Molecular weight 692.52 g/mol
- ✓ Molecular formula  $[\text{Fe}(\text{C}_{12}\text{H}_8\text{N}_2)_3]\text{SO}_4$
- ✓ Density 1.006 g/cm<sup>3</sup>
- ✓ Flash point does not flash
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 50mL
  
- ✓ H Statements H412
- ✓ P Statements P273

## Sales unit :

50mL Ferrouin indicator solution

---