

Standard Analysis Certificate

PRODUCT :n-PENTANE 95 % PUREX ,analytical grade
CODE :P06405

TEST	U.M.	SPECIFICATION
Refractive index at 20°C	-	1.355 - 1.359
Water content (K.F.)	mg/Kg	<= 150
Non volatile residue	mg/Kg	<= 10
Colour	Hazen	<= 10
Aromatic compounds	mg/Kg	<= 20
Assay (GC)	%	>= 95
n-hexane	%	<= 0.4
Cyclopentane	%	<= 1
2,2-dimethylbutane	%	<= 1
Total sulphur (S)	ppm	<= 2
Free acid (as CH ₃ COOH)	mg/Kg	<= 10

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE
P. GUERARD (PEYPIN)