



**HANDZENTRIFUGE**  
Manual Centrifuge

**Hettich**

ZENTRIFUGEN

# HANDZENTRIFUGE



## Maintenance-free, bench-mounted manually operated centrifuge

This convenient centrifuge can be easily clamped to a bench or table, and operated without electricity. Its self-lubricating gear makes it virtually maintenance-free.

It can be operated with two different rotors.

With rotor **1014** conical 15 ml tubes (Cat. No. 0515, non-graduated or No. 0516, graduated) can be accelerated to 1,298 RCF.

Rotor **1025** also accommodates a variety of other commercial tubes (also with round bottom) up to 15 ml. The tubes may have a max. length of 125 mm and a max. diameter of 17 mm.



illustrated with rotor 1014  
↘90°

illustrated with rotor 1025  
↘50°

ROTORS	4-place	4-place
Max. capacity	4 x 15 ml <sup>1)</sup>	4 x 15 ml
Max. RPM (speed)	3,000 min <sup>-1</sup>	3,000 min <sup>-1</sup>
Max. RCF	1,298	1,077
<b>Cat. No.</b>	<b>1014</b>	<b>1025</b>

<sup>1)</sup> conical only

HANDZENTRIFUGE	
<b>Manual centrifuge, without rotor</b>	
Dimensions (HxWxD)	285 x 140 x 175 mm
Weight	approx. 0.8 kg
<b>Cat. No.</b>	<b>1011</b>

LAB TECHNOLOGY

### Andreas Hettich GmbH & Co. KG

Föhrenstr. 12  
D-78532 Tuttlingen  
Germany  
www.hettichlab.com  
info@hettichlab.com  
service@hettichlab.com

Phone +49 (0)7461/705 -0  
Fax +49 (0)7461/705 -1122

**National Sales:** -12 00  
**International Sales:** -12 01  
**National Service:** -12 02  
**International Service:** -12 03