



HANDZENTRIFUGE
Manual Centrifuge

Hettich

ZENTRIFUGEN

HANDZENTRIFUGE



Maintenance-free, bench-mounted manually operated centrifuge

This convenient centrifuge can be easily clamped to a bench or table, and operated without electricity. Its self-lubricating gear makes it virtually maintenance-free.

It can be operated with two different rotors.

With rotor **1014** conical 15 ml tubes (Cat. No. 0515, non-graduated or No. 0516, graduated) can be accelerated to 1,298 RCF.

Rotor **1025** also accommodates a variety of other commercial tubes (also with round bottom) up to 15 ml. The tubes may have a max. length of 125 mm and a max. diameter of 17 mm.



illustrated with rotor 1014
↵90°

illustrated with rotor 1025
↵50°

ROTORS	4-place	4-place
Max. capacity	4 x 15 ml ¹⁾	4 x 15 ml
Max. RPM (speed)	3,000 min ⁻¹	3,000 min ⁻¹
Max. RCF	1,298	1,077
Cat. No.	1014	1025

¹⁾ conical only

HANDZENTRIFUGE	
Manual centrifuge, without rotor	
Dimensions (HxWxD)	285 x 140 x 175 mm
Weight	approx. 0.8 kg
Cat. No.	1011

LAB TECHNOLOGY

Andreas Hettich GmbH & Co. KG

Föhrenstr. 12
D-78532 Tuttlingen
Germany
www.hettichlab.com
info@hettichlab.com
service@hettichlab.com

Phone +49 (0)7461/705 -0
Fax +49 (0)7461/705 -1122

National Sales: -12 00
International Sales: -12 01
National Service: -12 02
International Service: -12 03