

Standard Analysis Certificate

PRODUCT :POTASSIO IODURO RPE-ACS
CODE :472730
METHOD :6600

TEST	U.M.	SPECIFICATION
Description	-	White or colourless solid
Identification	-	Positive
pH sol. 5% at 25° C	-	6.0 ÷ 9.2
Loss on drying at 150°C	%	<= 0.2
Water-insoluble matter	ppm	<= 50
Chloride + bromide (Cl)	ppm	<= 100
Iodate	ppm	<= 3
Phosphate	ppm	<= 10
Sulphate	ppm	<= 50
Heavy metals (Pb)	ppm	<= 5
Ba	ppm	<= 20
Ca	ppm	<= 20
Fe	ppm	<= 3
Mg	ppm	<= 10
Na	ppm	<= 50
Assay (oxidimetric)	%	>= 99.0

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)