

Standard Analysis Certificate

PRODUCT :POTASSIUM PERMANGANATE 0.1 N (0.02 mol/L) RPE**CODE** :473560000

TEST	U.M.	SPECIFICATION
Description Assay (potentiometry)	- N	Purple clear liquid 0.0998 - 0.1002

This volumetric solution is checked by sodium thiosulfate solution, tested with potassium dichromate volumetric standard traceable to the SRM 136e of NIST (National Institute of Standards and Technology)

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE**B. COULANGE (VDR)**