

# Standard Analysis Certificate

**PRODUCT** :POTASSIUM SODIUM TARTRATE RPE-ACS-ISO -Reag. Ph.Eur.- Reag. USP  
**CODE** :474110  
**METHOD** :9268

TEST	U.M.	SPECIFICATION
Description	-	White semitransparent crystals
Identification	-	Positive
pH sol. 5% at 25° C	-	6.0 ÷ 8.5
Ammonium	ppm	<=20
Chloride	ppm	<=10
Phosphate	ppm	<=20
Water-insoluble matter	ppm	<=50
Heavy metals (Pb)	ppm	<=5
Sulphate	ppm	<=50
Ca	ppm	<=50
Fe	ppm	<=10
Assay (non-aqueous medium)	%	99.0 ÷ 102.0

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)