

Standard Analysis Certificate

PRODUCT :L(-)PROLINE RPE
CODE :474700
METHOD :9280

TEST	U.M.	SPECIFICATION
Description	-	White crystalline powder
Identification	-	Positive
Potere rotat. spec. a 20°C (C=4; H2O)	°	-86.8 ÷ -84.5
pH 2,5% at 25° C	-	5.5 ÷ 7
Loss on drying	%	<= 0.3
Assay (non-aqueous medium)	%	99.0 ÷ 101.0

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE**B. COULANGE (VDR)**