

**40.401 CLASSIC HYGIENE PP GOWN**  
**40.705 2 pockets + snap fastener**
**FEATURES & BENEFITS**

- Single use hygiene gowns protect the working environment from particles coming from the body of the operator (dander, hair, sweating ...)
- They are suitable to HACCP procedures.
- They allow the operator to protect himself from the risks of soiling related to the activity he carries out

**TECHNICAL DESCRIPTION**

**Product name:** Classic hygiene PP gown  
**Product type:** Single use  
**Material:** Polypropylene  
**Confection :** Sewn and elasticated (collar ankles, wrists),  
**Snap fastener with pockets**

**Grammage :** 40 g/m<sup>2</sup>  
**Unit weight (L):** 72 g  
**Origin :** China

**40.401**


White


**40.705**


Dark blue



Size		Length (cm)	Width size (cm)
M	2	104	130
L	3	104	135
XL	4	104	140
XXL	5	104	145

**INDUSTRIES**

Industry / Chemical / Petrochemical / Food / Away from home food service / Non residential catering  
 Community / Administration / Sanitary / Maintenance / Cleaning

# 40.401 CLASSIC HYGIENE PP GOWN

## 40.705 2 pockets + snap fastener

### CERTIFICATIONS & STANDARDS

Manufactured under an ISO 9001 certified system

### STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

### LOGISTICAL INFORMATION



Outer case characteristics

Item	Size mm	Gross weight kg
40.401/02	470 x 280 x 280	4,13
40.401/03	470 x 280 x 280	4,64
40.401/04	470 x 280 x 280	4,93
40.401/05	470 x 280 x 280	5,38
40.401/06	470 x 280 x 280	5,8
40.705/03	470 x 280 x 280	4,24
40.705/04	470 x 280 x 280	4,3
40.705/05	470 x 280 x 280	4,53
40.705/06	470 x 280 x 280	4,54

Polybag characteristics

Size mm	Packaging
350 x 260	1 x 50 sachets individuels
350 x 260	1 x 50 sachets individuels
350 x 260	1 x 50 sachets individuels
350 x 260	1 x 50 sachets individuels
350 x 260	1 x 50 sachets individuels
350 x 260	1 x 50 sachets individuels
350 x 260	1 x 50 sachets individuels
350 x 260	1 x 50 sachets individuels
350 x 260	1 x 50 sachets individuels



EAN Codes

Item	Outer case	Polybag
40.401/02	3 700 381 505 596	3 700 381 505 602
40.401/03	3 700 381 500 829	3 700 381 502 274
40.401/04	3 700 381 500 836	3 700 381 502 281
40.401/05	3 700 381 500 843	3 700 381 502 298
40.401/06	3 700 381 506 722	3 700 381 506 739

Item	Outer case	Polybag
40.705/03	3 700 381 505 619	3 700 381 505 626
40.705/04	3 700 381 505 633	3 700 381 505 640
40.705/05	3 700 381 505 657	3 700 381 505 664
40.705/06	3 700 381 509 358	3 700 381 509 365