

40.301
40.301V
40.301R
40.701
CLASSIC HYGIENE PP GOWN WITH COLLAR
Snap fastener
FEATURES & BENEFITS

- Single use hygiene gowns protect the working environment from particles coming from the body of the operator (dander, hair, sweating ...)
- They are suitable to HACCP procedures.
- They allow the operator to protect himself from the risks of soiling related to the activity he carries out.

TECHNICAL DESCRIPTION

Product name: Classic hygiene PP gown with collar
Product type: Single use
Material : Polypropylene
Confection : Sewn and elasticated (collar ankles, wrists), zipper
Snap fastener
Grammage : 40 g/m²
Unit weight (M): 108 g
Origin : China

40.301V


Green


40.701


Dark blue

40.301


White


40.301R


Red



Size		Length (cm)	Width size (cm)
M	2	100	130
L	3	104	135
XL	4	104	140
XXL	5	104	145
XXXL	6	104	150


INDUSTRIES

Industry / Chemical / Petrochemical / Food / Away from home food service / Non residential catering / Community / Administration / Sanitary / Maintenance / Cleaning

CERTIFICATIONS & STANDARDS

Manufactured under an ISO 9001 certified system

40.301**40.301V****CLASSIC HYGIENE PP GOWN WITH COLLAR****40.301R****Snap fastener****40.701****STORAGE CONDITIONS**

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

LOGISTICAL INFORMATION

Outer case characteristics

Item	Size mm	Gross weight kg
40.301/02	470 x 280 x 280	4,07
40.301/03	470 x 280 x 280	4,5
40.301/04	470 x 280 x 280	4,9
40.301/05	470 x 280 x 280	5,3
40.301/06	470 x 280 x 280	5,8
40.301V/04	470 x 280 x 280	4,4
40.301V/05	470 x 280 x 280	4,8
40.301R/04	470 x 280 x 280	4,5
40.301R/05	470 x 280 x 280	4,4
40.301R/06	470 x 280 x 280	4,9
40.701/02	470 x 280 x 280	4,4
40.701/03	470 x 280 x 280	4,62
40.701/04	470 x 280 x 280	4,8
40.701/05	470 x 280 x 280	5,2

Polybag characteristics

Size mm	Packaging
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags
350 x 260	1 x 50 individual polybags

Codes EAN

Item	Outer case	Polybag	Item	Outer case	Polybag
40.301/02	3 700 381 505 510	3 700 381 505 527	40.301R/04	3 700 381 506 548	3 700 381 506 562
40.301/03	3 700 381 500 799	3 700 381 502 243	40.301R/05	3 700 381 506 555	3 700 381 506 579
40.301/04	3 700 381 500 805	3 700 381 502 250	40.301R/06	3 700 381 508 313	3 700 381 508 320
40.301/05	3 700 381 500 812	3 700 381 502 267	40.701/02	3 700 381 505 534	3 700 381 505 541
40.301/06	3 700 381 506 999	3 700 381 507 002	40.701/03	3 700 381 500 911	3 700 381 502 366
40.301V/04	3 700 381 505 558	3 700 381 505 565	40.701/04	3 700 381 500 928	3 700 381 502 373
40.301V/05	3 700 381 505 572	3 700 381 505 589	40.701/05	3 700 381 500 935	3 700 381 502 380

KOLMI HOPEN SASU

Bld de la chanterie | 49124 Saint Barthélemy d'Anjou – France

Bld de la chanterie | BP 10059 – 49181 Saint Barthélemy d'Anjou Cedex – France

T : 33 (0) 241 963 434 | F : + 33 (0) 241 963 453 | sales@kolmi-hopen.com | www.h-open.com