

40.550**HYGIENE PP+PE LAB COAT 60g/m²****MICROPOROUS ANTISTATIC WITH COLLAR + SNAP FASTENER****FEATURES & BENEFITS**

- Pleasant and lightweight, polypropylene is lined with a PE film for a better a better seal.
- The PE film binds PP particles and reduces particle release from fibers.
- Model with flapless collar and elasticated cuffs.
- Suitable for protected environments.
- Individually packaged for controlled hygiene and distribution.
- Pocketless to prevent objects from falling out.
- Tested in accordance with EN 1149-5 : 2018, but does not qualify as personal protective equipment.

**TECHNICAL DESCRIPTION**

Product name:	Hygiene Lab Coat Microporous Antistatic
Product type :	Single use, non sterile
Material :	PP+PE
Confection :	Sewn & elasticated (collar & wrists) Snap fastener
Color :	White
Grammage :	60 g/m ² (±10%)
Origin :	China

Item	Size	Length (± 2 cm)	Width (± 2 cm)
40.550/02	M	104 cm	125 cm
40.550/03	L	104 cm	135 cm
40.550/04	XL	104 cm	140 cm
40.550/05	XXL	104 cm	145 cm

INDUSTRIES

Industrial / Restaurant & catering / Community services /
Hygiene

CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001.

40.550

HYGIENE PP+PE LAB COAT 60g/m²

MICROPOROUS ANTISTATIC WITH COLLAR + SNAP FASTENER

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTIC INFORMATIONS

OUTERCASE SPECIFICATIONS					CDT SPECIFICATIONS	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	CDT	EAN CODE
40.550/02	500 x 300 x 320	6,9	36 outercases	3 700 381 509 341	1 x 50 Individual polybags	3 700 381 509 334
40.550/03		7,3		7 400 763 014 001		7 400 763 014 032
40.550/04		7,7		7 400 763 014 046		7 400 763 014 070
40.550/05		8,1		7 400 763 014 093		7 400 763 014 107