

<b>TECHNICAL DATA SHEET</b>	
<b>OS1-1</b>	

<b>PRODUCT</b>
----------------

Inoculating loop, 1  $\mu$ m  
Sterile R

<b>CHARACTERISTICS</b>
------------------------

RAW MATERIAL	Natural Polystyrene
DIMENSIONS	length 194 mm, weight 1 g
NOMINAL CAPACITY	1 $\mu$ m

<b>PRODUCT QUALITY</b>
------------------------

PRODUCTION SYSTEM	Controlled production process aseptic Electron Beam Sterilization
SHELF LIFE	3 years from date of Production

<b>PACKING</b>
----------------

SETTING	1 piece in a sleeve, 250 pcs. per bag, 2000 pcs. per carton
STORING CONDITIONS	room temperature

<b>CARTON DIMENSION</b>
-------------------------

* LENGTH	600 mm
* WIDTH	400 mm
* HEIGHT	290 mm