

TECHNICAL DATA SHEET OS10-1
--

PRODUCT

Inoculating loop, 10 μ m
Sterile R

CHARACTERISTICS

RAW MATERIAL	Natural Polystyrene
DIMENSIONS	length 195 mm, weight 1 g
NOMINAL CAPACITY	10 μ m

PRODUCT QUALITY

PRODUCTION SYSTEM	Controlled production process aseptic Electron Beam Sterilization
SHELF LIFE	3 years from date of Production

PACKING

SETTING	1 piece in a sleeve, 250 pcs. per bag, 2000 pcs. per carton
STORING CONDITIONS	room temperature

CARTON DIMENSION

* LENGTH	600 mm
* WIDTH	400 mm
* HEIGHT	290 mm