

TECHNICAL DATA SHEET

PC122PPS

Conical container 125 ml
Twisted cup

CHARACTERISTICS

* RAW MATERIAL	Natural Polypropylene
* DIMENSIONS	∅60mm H 75mm /
* CAPACITY	125 ml

PACKING

* SIZE OF BAG	1000x600 mm
* SETTING	200 pcs in a bag, 200 pcs
* STORING CONDITIONS	Room temperature

PRODUCT QUALITY

* PRODUCTION SYSTEM	Controlled production process ASEPTIC
* SHELF LIFE	3 years from date of Production

CASE DIMENSION

* LENGTH	600 mm
* WIDTH	400 mm
* HEIGHT	290 mm

NOEX Spółka z ograniczoną odpowiedzialnością Sp.j
Fabianowska 128, Komorniki