

<b>Product :</b>	Steam Autoclave Sterilization Process Indicator Tape REF: H-ATDA	
<b>Description :</b>	This tape is designed as a process indicator for steam autoclaving (steam sterilization). Recommended for use with disposable fabrics (mainly non-woven fabrics). Black lines will indicate required steam sterilization conditions have been attained.	
	<u>Typical Value</u>	<u>Range</u>
<b>Thickness :</b>	0.165mm (6.5 mils)	0.155 – 0.175mm
<b>Adhesion :</b>	8 N / 25mm (30 oz/in)	6-10 N / 25mm
<b>Tensile Strength :</b>	90 N / 25mm (20 lbs/in)	75 – 95 N / 25mm
<b>Elongation :</b>	7%	6 – 8%
<b>Temperature range :</b>	As required for sterilization	
<b>Meets Specification No.</b>	ANSI/AMMI ST-60 (Class 1), BSI-7893, EN ISO 11140-1 Class 1 C.I.D. A A-54589A (US)	
<p>The above specification represents average values. Actual production can vary. The selection and application of a tape is the responsibility of the buyer. Tape should be stored at room temperature and should not be exposed to sunlight. We reserve the right to make changes, without notice, to our products and specifications.</p>		
May 2010		