

EU Declaration of Conformity

Product	Laboratory Equipment	File Number	TFR909081_001
Manufacturer	Cole-Parmer Ltd Beacon Road Stone, Staffordshire ST15 0SA United Kingdom		TFR909081_002

This declaration of conformity is issued under the sole responsibility of the manufacturer

Object of Declaration Shaker Incubator
(reference the attached list of catalogue numbers)

The object of the declaration described above is in conformity with the relevant Union Harmonisation Legislation:

Low Voltage Directive	2014/35/EU
EMC Directive	2004/108/EC
RoHS Directive	2011/65/EC

References to the relevant harmonised standards used or references to the other technical specifications in relation to which conformity is declared:

IEC/EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1: General requirements.
IEC/EN 61010-2-010:2003	Particular requirements for laboratory equipment for the heating of materials.
IEC/EN 61010-2-051:2015	Particular requirements for laboratory equipment for mixing and stirring.
IEC/EN 61326-1:2006	Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1: General requirements (Class A).

Signed for and on behalf of the above manufacturer

Additional Information	Year of CE Marking: 2013
Place of Issue	Stone, Staffordshire, UK
Date of Issue	29 April 2013. Revised 09 March 2017
Authorised Representative	Steve Marriott
Title	Technical Director
Signature	

Catalogue Number	Model	Description
79520-00	SI-200D-120	Shaking Incubator, 51 L capacity, 120 V, 60 Hz
79520-04	SI-200D	Shaking Incubator, 51 L capacity, 230 V, 50 Hz
79520-08	SIM-200D-M-120	Microplate Shaking Incubator, 120 V, 60 Hz
79520-10	SIM-200D-M	Microplate Shaking Incubator, 230 V, 50 Hz
79520-14	SI-200D-L-120	Large Capacity Shaking Incubator, 115 L capacity, 120 V, 60 Hz
79520-16	SI-200D-L	Large Capacity Shaking Incubator, 115 L capacity, 230 V, 50 Hz
79520-20	SIC-200D-C-120	Large Capacity Refrigerated Shaking Incubator, 115 L capacity, 120 V, 60 Hz
79520-22	SIC-200D-C	Large Capacity Refrigerated Shaking Incubator, 115 L capacity, 230 V, 50 Hz