

Rack System Combifix SM A



Description

Basic rack with 5 clamping strips h. For fastening different upright vessels. With additional strips also suitable for other vessels.

Details

Order No.: 0050 154

Accessory: Clamping Strip h 0050 400

Clamping strip v 0050 399

Cramp strip 0050 401

Spring strip 0050 402

Test tube rack with 68 drillings à 14 mm Ø

Test tube rack with 44 drillings à 16 mm Ø

Test tube rack with 44 drillings à 18 mm Ø

Test tube rack with 14 drillings à 30 mm Ø