

Universal tray SM



Description

Coated tray with drillings for individual fastening of spring clamps, test tube racks, or holders for 1 microplate. Shaking platform 560 x 400 mm.

Loading Capacity Universal Tray SM Spring clamps (stainless steel)

Size (ml)	10	25	50	100	250	500	1000	2000	3000	5000
Pieces (max.)	117	64	63	49	23	15	11	6	4	2

Loading Capacity Universal Tray SM Test tube racks (with hinged foot)

Pieces max: 5

Please note:

For fastening test tube racks on a universal tray **test tube rack** and **hinged foot** are necessary.

Details

Order No.: 0051 472

Accessory: Spring clamps 10 ml
Spring clamps 25 ml
Spring clamps 50 ml
Spring clamps 100 ml
Spring clamps 250 ml
Spring clamps 500 ml
Spring clamps 1000 ml
Spring clamps 2000 ml
Spring clamps 3000 ml
Spring clamps 5000 ml
Test tube rack with 68 drillings à 14 mm Ø
Test tube rack with 44 drillings à 16 mm Ø
Test tube rack with 44 drillings à 18 mm Ø
Test tube rack with 14 drillings à 30 mm Ø
Hinged foot for test tube racks