

# TECHNICAL DATASHEET

## 3580010 Puller stainless steel 25 cm

Other tools for clean rooms

**Individualization on request!**

### Product description

Autoclavable puller made of stainless steel.

- o For glass cleaning
- o Autoclavable



**MADE IN EUROPE**

<b>ISO</b>	3	4	5	6	7	8	9
<b>GMP</b>			A/B		C	D	

### Product specifications

Article name | Article number  
Bemaßung | Weight  
Material

Puller stainless steel 25 cm | 3580010

25 x 17 x 3 cm | 192 g

- o Halter: Edelstahl 1.4016
- o Griff: Gummi
- o Gummilippe: Kautschuk

Packing unit | Selling unit

1 Pieces | 1 PU

Quantity in carton

10 Pieces

Carton size | Carton weight

58 x 16 x 5 cm | 0,3 kg

Pallet

120 x 80 x 220 cm

Customs tariff number

87169090

Country of origin

Italy

EAN-Number

4065237022771

Compatible accessories

Our sales department will be happy to advise you on suitable accessories

### Industry Sectors



This classification is a non-binding recommendation. The user has to approve the individual suitability and compatibility for his individual application area. The condition of the material is to check regularly, damaged products have to be exchanged promptly. All information without warranty. Subject to changes and deviations.