

TECHNICAL DATA SHEET

Part Number	501198
Description	Measuring pipette with attached syringe



Benefits :

- ✓ With attached syringe for simple and risk-free pipetting
- ✓ With durable brown graduation
- ✓ Zero point at the top, graduated to the tip

Certificate :

- ✓ Non sterile

Technical specifications :

- ✓ Contents 10 ml
- ✓ Graduation 0,1 ml
- ✓ Outer diameter 1 mm
- ✓ Length 385 mm

Sales unit :

1 Measuring pipette with attached syringe