

**Bellco Glass, Inc.**  
**Specification Sheet**  
**Product #: 1965-11500**

**Description: Microcarrier Spinner Flask Only**

**Size: 15000mL**

**Material: Type 1 Class A Borosilicate Glass**

**Dimensions: 294mm OD x 460mm Overall Height, 358mm across sidearms**

**Two 45mm Sidearms angled 30° offset center**


**100mm Screw Thread Center Neck**

**Flat Bottom – no dimple**

**Graduations: 3L to 15L, every 1L – White Ink (Approx Volumes)**

**Supply with: Two 45mm PP Blue Screw Caps & PP Blue Pouring Rings**

**Application / Design: Micro Carrier Spinner Flask is designed for the efficient growth of suspension cultures and growth of cells on Microcarrier beads**



**The information contained in this specification sheet we believe is correct to the best of our knowledge. The recommendations and suggestions herein are made without guarantee or representation as to results. We recommend that adequate tests be made in your laboratory or plant to determine if this product meets all of your requirements.**



**BELLCO GLASS, INC.**



Bellco Glass, Inc. · 340 Edrudo Road · Vineland, NJ · 08360

Ver: 7/05