

Standard Analysis Certificate

PRODUCT :HEPTANE mixture of isomers PUREX ,analytical grade

CODE :524380

TEST	U.M.	SPECIFICATION
Appearance	-	Clear colourless liquid
Identification	-	Conform
Colour	Apha	<= 10
Refractive index at 20°C	-	1.387 - 1.407
Density at 15°C	g/l	690 - 725
Distillation range	°C	85 - 105
End/initial boiling points difference	°C	<= 10
Water content (K.F.)	mg/Kg	<= 100
Non volatile residue	mg/Kg	<= 10
Aromatic compounds	mg/Kg	<= 100
n-hexane	%	None
Toluene	mg/Kg	None

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)