

# Standard Analysis Certificate

**PRODUCT** :PHOSPHATE BUFFER pH 7.4  
**CODE** :PS0740

TEST	U.M.	SPECIFICATION
pH	unité pH	7.35 - 7.45
Temperature of measurement	°C	19 - 21

Accuracy : +/- 0.05 pH units - Standardized against NIST reference materials.

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

**QUALITY CONTROL RESPONSIBLE**  
**P. GUERARD (PEYPIN)**