



KENTRON
MICROBIOLOGY

"The A-Box, a square anaerobic jar, that makes your anaerobic cultures much simpler."



*Looking for a
square anaerobic jar?
Try the A-Box*



KENTRON
MICROBIOLOGY

*"Make your lab-life easier
with the A-box."*

Looking for a square anaerobic jar?

If you are looking for a square anaerobic jar that fits microtiter plates for your microbiology cultures under anaerobic conditions, then look no further.

We proudly present to you the A-Box. This unique product is designed to use with microtiter plates (deep wells) for your cultures under anaerobic conditions, large testtubes or petri-dishes. The A-box can be used with multiple anaerobic systems, Anaerobic Chambers and GasPak methods.

Finally the moment we have all been waiting for, it's here. The A-Box, a square anaerobic jar, that makes your anaerobic cultures much simpler.



Technical specifications:

External dimensions: 18.0 x 18.4 x 18.6 cm

Internal dimensions: 14.0 x 17.0 x 13.1 cm

Housing: PMMA Lid: stainless steel, lid contains 2 clips for catalysts. The A-box will be supplied with 1 rack for microtiterplates (test tube racks optional).



Kentron Microbiology B.V.
Logistiekweg 20-27
7007 CJ Doetinchem The Netherlands
Tel: +31-314-380527
Fax: +31-314-382899
Email: office@kentron-microbiology.nl
Website: www.kentron-microbiology.nl