

## Data sheet

**Description:** Volumetric flask 5 ml with poly stopper NS 10/ 19

**Code no.:** 296 01 53

---

Material:

DURAN (Borosilicat glass 3.3)

Chemical properties:

- Hydrolytic resistance at 98°C: Class 1 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 1 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Physical properties:

- Permissible operation temperature: max. 500°C
- Transformation temperature: 525°C
- Density: 2,23 g/cm<sup>3</sup>
- Linear expansion coefficient: 3,3 x 10<sup>-6</sup> K<sup>-1</sup>

**Height:**

- 70 +- 5 mm

**Tolerance:**

- +- 0,04 ml
- 

Hirschmann Laborgeräte  
GmbH & Co.KG



Roland Benz  
Quality Manager

Eberstadt, 04.09.2023