

TECHNICAL DATA SHEET

Part Number	501163
Description	Measuring pipette



Benefits :

- ✓ Measuring pipette, class B
- ✓ AR®-glass
- ✓ Brown graduation with silver diffusion ink durable
- ✓ Zero point at top
- ✓ Graduated to tip
- ✓ Calibrated for EX

Certificate :

- ✓ DIN EN ISO 835

Technical specifications :

- ✓ Contents : 1 ml
- ✓ Length +/- 5 mm : 360 mm
- ✓ Graduation : 0,01 ml
- ✓ Tolerance : 0,01

Sales unit :

12 Measuring pipette
