

TECHNICAL DATA SHEET

Part Number	501165
Description	Measuring pipette



Benefits :

- ✓ Made from AR-glass
- ✓ With durable brown graduation
- ✓ Zero point at the top

Certificate :

- ✓ Class B
- ✓ DIN EN ISO 835

Technical specifications :

- ✓ Contents 2 ml
- ✓ Graduation 0,01 ml
- ✓ Tolerance 0,015
- ✓ Length 360 mm (± 5 mm)

Sales unit :

12 Measuring pipettes