

Technical data sheet

Article: Graduated pipette 5:0,05 ml class B

Code no.: 100 01 54

Material:

AR-Glas (Soda-lime glass)

Chemical properties:

- Hydrolytic resistance at 98°C: Class 3 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 3 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Tolerance:

± 0,045 ml

Standard:

ISO 835

Hirschmann Laborgeräte
GmbH & Co.KG



Roland Benz
Quality Manager

Eberstadt, 17.06.2019