

## TECHNICAL DATA SHEET

<b>Part Number</b>	501170
<b>Description</b>	Measuring pipette



### Benefits :

- ✓ Made from AR-glass
- ✓ With durable brown graduation
- ✓ Zero point at the top

### Certificate :

- ✓ Class B
- ✓ DIN EN ISO 835

### Technical specifications :

- ✓ Contents 10 ml
- ✓ Graduation 0,1 ml
- ✓ Tolerance 0,075
- ✓ Length 360 mm ( $\pm 5$  mm)

### Sales unit :

12 Measuring pipettes