

Bellco Glass, Inc.
Specification Sheet
Product #: 1965-20001R

Description: SST 2-Port Cap Assembly

Size: 1L

For Bellco Microcarrier Spinner Flask with angled 45mm sidearms

Material: 316L SST with Electropolish Finish of Ra32 μ or better – inside & out

Tubing Specifications: .250" OD x .152" ID


Tube 1: 1.0" long above & below cap slug

Tube 2: 1.0" long above cap slug & angled to extend down close to bottom as possible of a Bellco Microcarrier 1L Spinner Flask

Seam weld – top & bottom side

Supplied with blue Polypropylene open top 45mm screw cap & white silicone gasket

Application/usage: Media Addition, cell incubation, sub-culture, bioreactor inoculation



The information contained in this specification sheet we believe is correct to the best of our knowledge. The recommendations and suggestions herein are made without guarantee or representation as to results. We recommend that adequate tests be made in your laboratory or plant to determine if this product meets all of your requirements.



BELLCO GLASS, INC.



Bellco Glass, Inc. · 340 Edrudo Road · Vineland, NJ · 08360

Ver: 7/05