

**Bellco Glass, Inc.**  
**Specification Sheet**  
**Product #: 1965-20003R**

**Description: SST 2-Port Cap Assembly**

**Size: 3L**

**For Bellco Microcarrier Spinner Flask with angled 45mm sidearms**

**Material: 316L SST with Electropolish Finish of Ra32 $\mu$  or better – inside & out**

**Tubing Specifications: .250" OD x .152" ID**


**Vent Port: 1.0" long above & below cap slug**

**Media Port: 1.0" long above cap slug & angled to extend down close to bottom as possible of a Bellco Microcarrier 3L Spinner Flask**

**Seam weld – top & bottom side**

**Supplied with blue Polypropylene open top 45mm screw cap & white silicone gasket**

**Application/usage: Media Addition, cell incubation, sub-culture, bioreactor inoculation**



**\*Blueprint revision 0 (6/26/2007)**

**The information contained in this specification sheet we believe is correct to the best of our knowledge. The recommendations and suggestions herein are made without guarantee or representation as to results. We recommend that adequate tests be made in your laboratory or plant to determine if this product meets all of your requirements.**



USA

BELLCO GLASS, INC.



Bellco Glass, Inc. · 340 Edrudo Road · Vineland, NJ · 08360

Ver: 7/05