

Bellco Glass, Inc.
Specification Sheet
For Product #: 1965-51100

Description: 1L Water Jacketed Micro Carrier Spinner Flask Only

Micro Carrier Spinner Flask is designed for the efficient growth of suspension cultures and growth of cells on microcarrier beads. The flasks incorporate a fluid-dynamic impeller assembly to assure optimal lift of suspension cells or microcarrier beads at very low mixing speeds without damage to cells or beads.

Flask Specifications:


Material: Type 1 Class A Borosilicate Glass

Dimensions: 170mm OD x 260mm Overall Height

Two 45mm Sidearms angled 60° offset center

One 100mm Screw Thread Center Neck

Supply with: One 100mm Cap, Two 45mm PP Blue Screw Caps & PP Blue Pouring Rings



The information contained in this specification sheet we believe is correct to the best of our knowledge. The recommendations and suggestions herein are made without guarantee or representation as to results. We recommend that adequate tests be made in your laboratory or plant to determine if this product meets all of your requirements.



BELLCO GLASS, INC.



Bellco Glass, Inc. · 340 Edrudo Road · Vineland, NJ · 08360

Ver: 7/05