


**Bellco Glass, Inc.**  
**Specification Sheet**  
**Product #: 2029-00018**

**Description: Black Phenolic Rubber Lined Screw Cap**  
**Size: 18-415 / Deep Skirted**  
**Material of Cap: Black Phenolic**  
**Dimensions of Cap: 0.807" OD x 0.70" Long, approx.**  
**Material of Liner: 14B White Rubber**  
**Liner is bonded to cap**  
**Autoclavable**  
**Designed for culture work**  
**For all Laboratory Bottles that have a 18mm GCM I Thread**



**The information contained in this specification sheet we believe is correct to the best of our knowledge. The recommendations and suggestions herein are made without guarantee or representation as to results. We recommend that adequate tests be made in your laboratory or plant to determine if this product meets all of your requirements.**



**BELLCO GLASS, INC.**



Bellco Glass, Inc. · 340 Edrudo Road · Vineland, NJ · 08360

Ver: 7/05