

Datenblatt - Data Sheet

Customer:

Article:

Laborflasche GL45 - Rasotherm

Media Bottle GL45 - Rasotherm

Article number:

95206001 - 95206007

Material:

Glas / glass: Boro 3.3 ISO 3585; autoclavable


Kappe / cap: Polypropylene; autoclavable; max. 140°C max

Dimension/Standard:

Kappe / cap: GL45 – Gewinde / GL45 screw thread

Glas / glass: ISO 4796

Article	Size	Overall height	Outside diameter
	[mL]	[mm]	[mm]
95206001	100	100	Ø56
95206002	250	138	Ø70
95206003	500	176	Ø86
95206004	1000	225	Ø101
95206005	2000	260	Ø136
95206006	5000	330	Ø181
95206007	10000	410	Ø227



Glass properties:

Water resistance tested to DIN ISO 719 (at 98 °C): HGB 1

Water resistance tested to DIN ISO 720 (at 121 °C): HGA 1

Alkali resistance tested to DIN ISO 695: class A2

Acid resistance tested to ISO 1776: $\text{Na}_2\text{O} \leq 100 \mu\text{g}/\text{dm}^2$