

Description

- Boro glass tube 3.3 straight edge round bottom 12x75mm thickness 0.9-1.0mm

Part number

- 511422

Pieces by sales : 100



Description

- Test tubes, Borosilicate glass 3.3
- Diam : 12mm
- Height: 75mm
- Wall thickness: 0.8 - 1.0mm
- With rim: without
- Catalogue Number 4.658 114
- Manufacturer : Scherf Prazision
- Manufacturer P/n 90106015
- Material : Boro 3.3
- With beaded rim : No
- Length mm 75
- Ø : 12.0 mm
- Wall thickness mm : 0.8 - 1.0
- Dimensions (ØxL) : 12 x 75 mm