

DECLARATION OF CONFORMITY

Cat.No. FG-NP01

Year of CE marking: 2020

Nippon Genetics Europe GmbH
Mariaweißerstraße 28-30
D-52349 Düren, Germany
Phone: +49 2421 554960
Fax: +49 2421 5549611

Product name: FastGene® Photometer
Model name: FG-NP01 FastGene® NanoSpec Photometer

This is to certify that the above listed product was tested with the required criteria levels of the follow-mentioned Generic Standards or Product Family Standards and/or Basic Standards based-on the essential conformity requirements of EMC Directive of 2014/30/EC and of Low Voltage Directive 2014/35/EC.

Standards: EN 61326-1:2013
EN 55011:2016/A1:2017
EN 61010 – 1:2010

Name and address of the person responsible for documentation:

Dr. Manuel Franke
Name

Mariaweißerstraße 28-30, 52349 Düren
Address

Identification of the signatory:

Dr. Jürgen Lünzer, CEO of Nippon Genetics Europe
Name, Function

02.03.2021 [Signature]
Place, Date, Signature