

Standard Analysis Certificate

PRODUCT :ACIDO SOLFORICO 96% RPE - ISO
CODE :410300
METHOD :6585

TEST	U.M.	SPECIFICATION
Description	-	Clear oily liquid
Colour (APHA)	-	<=10
Colour of 2N solution(APHA)	-	<=10
Identification	-	Positive
Density at 20° C	-	1.834 ÷ 1.836
Ammonium	ppm	<=1
Chloride	ppm	<=0.1
Heavy metals (Pb)	ppm	<=0.8
Nitrate	ppm	<=0.2
Residue on ignition	ppm	<=4
Subst. reducing KMnO4	ppm	<=2
As	ppm	<=0.01
Ca	ppm	<=0.2
Cd	ppm	<=0.005
Cu	ppm	<=0.01
Fe	ppm	<=0.1
K	ppm	<=0.1
Li	ppm	<=0.02
Mg	ppm	<=0.2
Na	ppm	<=0.5
Ni	ppm	<=0.05
Pb	ppm	<=0.02
Se	ppm	<=3
Sr	ppm	<=0.02
Zn	ppm	<=0.05
Assay (acidimetric)	%	95 ÷ 97

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)