

Standard Analysis Certificate

PRODUCT :METHANOL RS For HPLC Isocratic Grade
CODE :525100
METHOD :0001

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Assay GLC	%	>= 99.9
Acidity (Formic Ac)	ppm	<=20
Density at 20°C	-	0.791 ÷ 0.793
Refractive index at 20°C	-	1.3278 ÷ 1.3298
Boiling point	°C	63.6 ÷ 65.6
Residue on evaporation	ppm	<= 10
Water (KF)	ppm	<= 500
Transmittance UV at 210 nm	%	>=20
Transmittance UV at 235 nm	%	>= 80
Transmittance UV at 260 nm	%	>= 98

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)