

# Technical Specification

**PRODUCT : Potassium dichromate 0,04 mol/l (0,24 N) in 80 g/l HgSO<sub>4</sub> RS - For environmental analysis (COD determination)**

**CODE : 526711**

TEST	U.M.	SPECIFICATION
Assay	mol/L	0.0398 - 0.0402

DOMINIQUE DUTSCHER SAS

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in