

# Standard Analysis Certificate

**PRODUCT** :TOLUENE RE PURO  
**CODE** :386000  
**METHOD** :7960

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20° C	-	0.862 ÷ 0.872
Refractive index at 20°C	-	1.491 ÷ 1.501
Boiling point	°C	109.9 ÷ 111.4
Residue on evaporation	ppm	<=50
Water (K.F.)	ppm	<=300
Benzene	%	<= 0.02
Total sulphur	ppm	<=100
Assay (GLC)	%	>=99.5
Ethyl benzene + xylene	ppm	<= 750

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)