

# Standard Analysis Certificate

**PRODUCT** :XYLENE MIX OF ISOMERS RE PURO LOW IN BENZENE  
**CODE** :528250  
**METHOD** :20045

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Identity (IR)	-	Positive
Density at 20°C	-	0.862 ÷ 0.872
Refractive index at 20°C	-	1.4917 ÷ 1.5017
Boiling point	°C	137.5 ÷ 139.5
Residue on evaporation	ppm	<= 100
Water (K.F.)	ppm	<= 200
Benzene	ppm	<= 50
Total sulphur	ppm	<= 100
Assay (isomeric mix)	%	>= 98.5

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE
**B. COULANGE (VDR)**