

Standard Analysis Certificate

PRODUCT :CHLOROFORME STAB AMYLENE RE PURE
CODE :528325
METHOD :20018

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Color	APHA	<10
Identity (IR)	-	Positive
Density at 20°C	-	1.478 ÷ 1.488
Refractive index at 20°C	-	1.4461 ÷ 1.4491
Boiling point	°C	61 ÷ 61.5
Acidity (HCl)	ppm	<10
Residue on evaporation	ppm	<=20
Water (K.F.)	ppm	<= 300
Assay (GLC)	%	>= 99.9
Amylene	ppm	<= 60

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)