

# Standard Analysis Certificate

**PRODUCT** :Ethanol 95° denaturated RE According to European regulation**CODE** :528770

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<= 10
Density at 20°C	-	0.804 ÷ 0.812
Assay (alcohol)	% (v/v)	94.8 ÷ 96.8

Regulation 2016/1867 of 20th October 2016

Denaturation by European process for 100L:

1 L of isopropyl alcohol  
1 L of methylethylketone  
1 g of Bitrex**Date**

:11/01/2020

**QUALITY CONTROL RESPONSIBLE****B. COULANGE (VDR)**