

Standard biopsy punches

The Standard biopsy punches make it possible to perform skin biopsies in various locations with great precision and minimal tissue trauma.

View



Benefits

Optimal results

The Biopsy Punch guarantees the highest quality due to the stainless steel blade with micro-grinding finish. It allows a representative cross-section with sufficient depth and width through the entire cutis.

User friendly

The size indication on each biopsy ensures an easy identification.

Ready to use

Each Punch is individually sterile packed and ready to use.

Ordering information

Standard biopsy punches

REF	Size	PU
48101	1,0 mm	20 pcs
48115	1,5 mm	20 pcs
48201	2,0 mm	20 pcs
48251	2,5 mm	20 pcs
48301	3,0 mm	20 pcs
48351	3,5 mm	20 pcs
48401	4,0 mm	20 pcs
48501	5,0 mm	20 pcs
48601	6,0 mm	20 pcs
48801	8,0 mm	20 pcs